

## Reaction Conditions - TexVet20

RXN#1				Cocktail Volume (uL)	1130
Reagent	[ ]	[ ] Units	Desired [ ]	Des'd [ ] Units	Amount to add
Buffer	10	X		1 X	113
dNTP's	2	mM Each		0.2 mM Each	113
MgCl2	50	mM		1.5 mM	33.9
TV17-F-6FAM	20	uM		0.3 uM	16.95
TV17-R	20	uM		0.3 uM	16.95
Taq	5	U/uL		0.05 U/uL	11.3
DNA	5	ng/uL		10 ng/rxn (2 uL/rxn)	226
Autoclaved ddH2O			to 1130 uL		598.9
					1130
Dilute by 3 ul of PCR product with 15 ul of water					
RXN#2				Cocktail Volume (uL)	1130
Reagent	[ ]	[ ] Units	Desired [ ]	Des'd [ ] Units	Amount to add
Buffer	10	X		1 X	113
dNTP's	2	mM Each		0.2 mM Each	113
MgCl2	50	mM		1.5 mM	33.9
RW417-F-6FAM	20	uM		0.4 uM	22.6
RW417-R	20	uM		0.4 uM	22.6
TV20-F-HEX	20	uM		0.07 uM	3.955
TV20-R	20	uM		0.07 uM	3.955
RW45-F-NED	20	uM		0.11 uM	6.215
RW45-R	20	uM		0.11 uM	6.215
RW212-F-HEX	20	uM		0.085 uM	4.8025
RW212-R	20	uM		0.085 uM	4.8025
Taq	5	U/uL		0.05 U/uL	11.3
DNA	5	ng/uL		10 ng/rxn (2 uL/rxn)	226
Autoclaved ddH2O			to 1130 uL		557.655
					1130
Don't Dilute Prior To De-salting					

Use PCR Cycles PROPDNR		
5 min	@	94
30 sec	@	94
1 min	@	49
1 min	@	72
45 min	@	60

Is a 10ul reaction volume

Use PCR Cycles PROPDNR		
5 min	@	94
30 sec	@	94
1 min	@	58
1 min	@	72
45 min	@	60

Is a 10ul reaction volume

## TextVet20 Scoring Criteria

- Score the +A, although it is difficult to tell because the +A is so preferentially amplified.
- Allele bins are even

